


Molequin

WATCH SIZING GUIDE

INSTRUCTIONS FOR DETERMINING LUG WIDTH

- 1 Print off this page at full size. Make sure that it's sized at 100%.
- 2 Ideally remove the bracelet or strap from your watch and then place over the watch diagram that correlates to yours.
- 3 The black box should fit perfectly between the lugs of your watch and determine lug width. So a 19mm measurement will require a 19mm strap.



INSTRUCTIONS FOR DETERMINING WRIST SIZE

- 1 Start by cutting along the lines marked with an A.
- 2 Make a small cut at B (just half way - don't cut all the way through).
- 3 When you're ready wrap the tape measure around your wrist and slide the end through the slit at B, which will help you read off your measurement.
- 4 Short Size - from 14.5cms to 17.5cm.
Standard Size - from 16.5cms to 19.0cms.

